Reactivos GPL

Store at: +2+8°C.

Barcelona, España

 ϵ

- Hemoglobin -

HEMOGLOBIN Drabkin, Colorimetric

Presentatión:

CONT: R 1 x 50 mL. Cod. SU020 Cod. SU021 CONT: R 4 x 5 mL.

Procedure

Hemoglobin g/dL = (A) Sample x 37.7

Quantitative determination of Hemoglobin.

Only for in vitro use in clinical laboratory (IVD)

TEST SUMMARY

oxidized by potassium Hemoglobin is ferricvanide metHemoglobin, which is converted into cyanometHemoglobin, by potassium cyanide.

The intensity of the color formed is proportional to the Hemoglobin concentration in the sample1

REAGENTS COMPOSITION

Hemoglobin 50x	Potassium ferricyanide Potassium cyanide Dihydrogen potassium phosphate	0.60 mmol/L 0.90 mmol/L 2 mmol/L
Hemoglobin CAL (Optional)	Hemoglobin Standard Animal origin	15 g/dL

PRECAUTIONS

HEMOGLOBIN Ref: SU020 / SU021

Harmful (Xn): R20/21/22: Harmful by inhalation, in contact with skin and if swallowed. S7: Keep container tingly close. S28.1: After contact with skin, wash immediately with plenty of water. S45: In case of accident or if you feel unwell, seek medical advice immediately.

Cyanide (poison): The amount of cyanide in the Reagent Concentrate (50x) is appreciably less than the minimum lethal dose for an adult. Gaseous hydrogen cyanide will be released on contact with acids.

REAGENT PREPARATION AND STABILITY

Working reagent (WR):

For 5 mL 4.9 mL of distilled water + 2 drops of Reagent

For 250 mL 245 mL of distilled water + 1 vial (5 mL) of Reagent Mix well.

Stability: 30 days at 2-8°C, protected from the sunlight.

All the components of the kit are stable until the expiration date on the label when stored at 2-8°C, protected from light and contamination prevented during their use.

Do not use reagents over the expiration date.

Signs of Reagent deterioration:

- Presence of particles and turbidity

- Blank absorbance (A) at $540 \text{ nm.} \ge 0.01$ All the reagents of the kit are stable up to the end of the indicated month and year of expiry. Store tightly closed at 2-8°C. Do not use reagents over the expiration date.

SPECIMEN

Venous or capillary blood¹

Use anticoagulants like EDTA, heparin or oxalate.

Stability of the sample: 1 week at 2-8°C.

MATERIAL REQUIRED BUT NOT PROVIDED

- Spectrophotometer or colorimeter measuring at 540 nm.
- Matched cuvettes 1.0 cm. light path.

General laboratory equipment.

TEST PROCEDURE

1.	Assay conditions:	
	Wavelength:	. 540 nm
	Cuvette:	light path
	Temperature	15-25°C

2 Adjust the instrument to zero with distilled water.

Pipette: 3

	Blank	Standard	Sample
WR (mL.)	5.0	5.0	5.0
Hemoglobin Cal. (μL.)		20	
Sample (μL.)			20

Mix and incubate for 3 min. at room temperature (15-25°C).

Read the absorbance (A) of the samples and calibrator, against the Blank. The colour is stable for at least 48 hours.

CALCULATIONS

With factor²:

With calibrator:

(A)Sample Hemoglobin g/dL = x 15 (Standard conc.) (A)Standard

QUALITY CONTROL

Each laboratory should establish its own Quality Control scheme and corrective actions.

REFERENCE VALUES

Men	14 - 18 g/dL ≅ 8.7 – 11.2 mmol/L
Women	12 - 16 g/dL ≅ 7.5 – 9.9 mmol/L

(These values are for orientation purpose).

It is suggested that each laboratory establish its own reference

CLINICAL SIGNIFICANCE

The Hemoglobin is a protein that contains iron and that the red color to the blood. The Hemoglobin is in red globules and it is the one in charge of oxygen transport by the blood from the lungs to weaves.

When the level of Hemoglobin appears underneath the normal levels is describing an anemia that can be of different origins: primary anemia,

cancer, pregnancy, renal diseases, and hemorrhages. If the Hemoglobin levels appear high it can be due to: cardiopathies, dehydratation and stays in places of much altitude^{1,5,6}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENT PERFORMANCE

 $\underline{\textit{Measuring Range}}\text{:}$ From detection limit 0f 0.1 g/dL. to linearity limit of 18 g/dL., under the described assay conditions.

If results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L. and multiply result by 2.

Precision:

	Intra-assay n= 20		Inter-ass	ay n= 20
Mean (g/dL)	8.00	15.2	7.81	15.1
SD	0.29	0.33	0.19	0.26
CV (%9	3.59	2.19	2.51	1.74

Sensitivity: 1 g/dL. = 0.027A

Accuracy:

Results obtained GPL reagents did not show systematic differences when compared with other commercial reagents.

The results of the performance characteristics depend on the analyzer used

INTERFERING SUBSTANCES

A list of drugs and other substances that could interfere has been reported by Young et. $a^{\rm 5.6}.$

NOTES

Use clean disposable pipette tips for its dispensation.

BIBLIOGRAPHY

- Franco R S. Hemoglobin. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1294-1296 and 418.
- 2. Van Kampen EJ et al. Standarization of Hemoglobinometry Clin. Chim 1961;6:
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995. 3.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed. AACC 2001 Burtis A. et al. Tietz Textbook of Clinical Chemistry, 3rd ed. AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed. AACC 1995.

